

## BINSON'S HOME HEALTH CARE CENTERS

**BINSON'S HOME HEALTH CARE CENTERS  
INCONTINENCE & MEDICAL SUPPLY PROGRAM**  
Monday through Friday from 8:00 A.M. to 5:30 P.M.  
• **Toll Free Number 1-888-217-9610**  
• TTY Relay Service at 711 For Hearing Impaired  
**24-HOUR EMERGENCY SERVICE AVAILABLE**

Welcome to Binson's State of Indiana Incontinence and Medical Supply Program. We are committed to providing you with quality incontinence products. All products will need to be ordered by your physician. We also offer an in-house nursing staff for telephone consultation and initial nursing assessments.

Beginning February 1, 2008, Binson's began providing bowel/bladder incontinence supplies by FedEx, United Parcel System (UPS) or United States Postal System (USPS) through a home delivery contract with OMPP. Your products will be shipped to you monthly unless you report changes.

Our team of customer service representatives is available to answer your questions and assist with your product needs. If you have questions regarding your supplies or did not receive your supplies, please call us at **1-888-217-9610**.

**PRODUCT SUPPLIED: DRAINAGE BAG, LEG BAG**

### **DESCRIPTION:**

- Available in latex and vinyl.
- Made of flexible, disposable material that allows for the collection and drainage of urine when using a catheter. Can be washed out and reused.

### **PRODUCT USE:**

#### **Drainage bag**

- Attach the bag to the bed side or chair side with the hook on the strap.
- Remove the cap from the drainage tube and connect it to the catheter without touching the catheter tip or tip of the tubing.
- Secure the tube to the bed sheet or chair side with the bed clip or strap.
- Drain when necessary into the toilet or other collection area by opening the clamped tube at the bottom of the bag. Re-clamp after emptying.

#### **Leg Bag**

- Insert the straps through the slots.
- Position the bag on the inside of the leg.
- Adjust the straps to fit securely.
- Drain when necessary into the toilet or other collection area by opening the clamped tube at the bottom of the bag. Re-clamp after emptying.

## BINSON'S HOME HEALTH CARE CENTERS

Page 2

Patient Education

Drainage bag, Leg Bag

### **NURSING/SKIN CARE:**

- Bag must be below the level of the bladder to drain. Check and empty the bag regularly to prevent overflow. Throw away bags and tubing.
- Call your physician if you notice any change in urine color, odor, any burning or pain when emptying your bladder.
- Antiseptic wipes or alcohol wipes can be used to clean connection tubing before attaching to the catheter and after disconnecting.

### **CLEANING:**

- Bags may be soaked with a solution of vinegar and water to decrease mineral deposit build-up.
- Mix four (4) ounces of vinegar with one (1) quart of water. Soak connectors and bag valves for six to eight hours, rinse with clear water.
- Repeat every three to four days.